

Technical Data

MULTIDOS Type 10004

Type of instrument	Field class therapy dosimeter according to DIN 6817 and IEC 731
Application	Dose and dose rate measurements (charge and current measurements) in radiation therapy for quality assurance, patient dosimetry and field analysis
Electrical safety class	IEC 601-1 and IEC 1010-1 (EN 60601-1 and EN 61010-1)
Number of channels	12
Sampling frequency	500 Hz
Warming up period	< 15 min
Internal adjustment time	< 1.5 s
Accuracy of the current and charge measurement	< ± 0.5%
Interval time	(6 ... 9999) s
Chamber voltage	+ 400 V ± 2%
Long term stability	< ± 0.5% per year
Linearity	< ± 0.5% according to IEC 731
Maximum integration time	18 h
Measuring range for dose rate measurements with 0.5% resolution (with semi-conductor probe type 9111)	Low: (0.4 · 10 ⁻³ ... 0.4) Gy/min High: (40 · 10 ⁻³ ... 40) Gy/min
Measuring range for dose rate measurements with 0.5% resolution (with ionization chamber type 31002)	Low: (0.14 ... 140) Gy/min
Measuring range for current measurements with 0.5% resolution	Low: (10 ⁻¹¹ ... 10 ⁻⁸) A High: (10 ⁻⁹ ... 10 ⁻⁶) A
Measuring quantities and units	Absorbed dose to water, air kerma [Gy] Absorbed dose rate to water, air kerma rate [Gy/min] exposure [R] exposure rate [R/min] Charge [C] Current [A]
Temperature range	(10 ... 40) °C
Range of air humidity	(10 ... 80)% rel. humidity; maximum 20 g/m ³
Range of air pressure	(700 ... 1060) hPa
Zeroing of the amplifier	automatically in approx. 30 s
Leakage current (zero drift)	< ± 5 · 10 ⁻¹⁴ A
Power supply	115 V / 230 V; (50 ... 60) Hz
Current limitation of the high voltage supply	0.5 mA by combined active and passive current limitation
Time constant for current and dose rate measurements	< 0.5 s
Weight	approx. 5.1 kg
Outer dimensions	Width: 25.9 cm; height: 11.5 cm; depth 32.6 cm

One multi channel dosimeter offers five applications:

- Dual channel dosimetry
- Afterloading dosimetry
- Total body dosimetry
- Constancy checks
- Dynamic field dosimetry



MULTIDOS

PTW-Freiburg is an accredited calibration laboratory of the German Calibration Service for dosimeters and sensitometers.

PTWFREIBURG

PHYSIKALISCH-TECHNISCHE WERKSTÄTTEN
DR. PYCHLAU GMBH

✉ PTW-FREIBURG
LÖRRACHER STRASSE 7
D-79115 FREIBURG, GERMANY
☎ +49-761-49055-0
FAX +49-761-49055-70
E-MAIL ptw@ptw.de
WEB <http://www.ptw.de>

Subject to change without notice 5/00 – 1B

PTWFREIBURG

► Overview

MULTIDOS is a multi channel dosimeter which combines the latest developments in microprocessor based electrometer technology with state-of-the-art electronics to produce a truly superior multi channel radiation therapy dosimeter. MULTIDOS is easily configured by the user for a wide variety and range of measuring tasks in radiation therapy which include:

- Dual channel dosimetry according to IEC 731
- Direct patient dose measurements during afterloading, external beam and total body treatments
- Quality assurance constancy checks including flatness, symmetry and energy stability
- Dosimetry and analysis of dynamic and irregular fields

MULTIDOS may be used as a dual, six or twelve channel dosimeter. With the optional MULTIDOS extender ME 48 you will have the ability of making 47+1 simultaneous measurements. 47+1 means 47 field measurements and 1 reference channel measurement for a total of 48 channels and 48 simultaneous measurements.

MULTIDOS can be used with ionization chambers and semi-conductor detectors. Chamber and detector calibration factors are easily stored in MULTIDOS and you can also enter correction factors, for example to correct for different beam qualities. If you are measuring with vented ion chambers, you can input temperature and air pressure and MULTIDOS will calculate and apply the air density correction factor.

MULTIDOS offers 5 different configurations and meets or exceeds the standards for field class electrometers as specified in DIN 6817 and IEC 731. With a sampling frequency of 500 Hz, you will find that all of the data you collect with MULTIDOS will always be accurate, reproducible and extremely reliable. You can depend on MULTIDOS when simultaneous measurements of multiple channels is a must.

MULTIDOS configurations

MULTIDOS applications	Measuring task
2-channel dosimeter	Constancy check of the energy by measuring the dose in different depths
6-channel dosimeter	Patient dosimetry during afterloading therapy (single probe for the bladder and fivefold probe for the rectum)
6-channel dosimeter	Constancy check of flatness, symmetry and energy, using a test object (solid state phantom)
12-channel dosimeter	Patient dosimetry with 12 detectors during total body irradiation
47+1-channel dosimeter (extender ME 48 required)	Analysis of dynamic fields with the optional linear array LA 48 and the MULTIDOS extender ME 48, using the therapy beam analysers PTW-MP3 or MP3-S

